

# EBP07046 - IMMUNIZING PEPTIDE

SKU: EBP07046



---

## SPECIFICATIONS

<b>Formulation</b>	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
<b>Unit Size</b>	100 µg
<b>Storage Instructions</b>	Shipped at ambient temperature, store at -20°C.
<b>Accession ID</b>	NP_001184151.1; NP_001184150.1; NP_006194.2; NP_001098101.1; NP_001184152.1; NP_001159371.1
<b>Immunogen</b>	Peptide with sequence QPEACVIDDRSPDT, from the C Terminus of the protein sequence according to NP_001184151.1; NP_001184150.1; NP_006194.2; NP_001098101.1; NP_001184152.1; NP_001159371.1.
<b>Product Comments</b>	The peptide was used in the production of EB07046 - Goat Anti-PDE4D Antibody. It is currently untested in blocking applications.
<b>Peptide Sequence</b>	QPEACVIDDRSPDT
<b>Shipping Instructions</b>	Ambient
<b>Product Grade</b>	<a href="https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png">https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png</a>

## DOCUMENTS

- [Data Sheet](#)