

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP07048 - IMMUNIZING PEPTIDE**

**SKU:** EBP07048

## **SPECIFICATIONS**

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$ 

**Unit Size** 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP\_005663.1 ID

Immunogen Peptide with sequence CVDDSQDYYVGKKN, from the internal region of the protein sequence according to NP\_005663.1.

The peptide was used in the production of EB07048 - Goat Anti-PSCA (aa69-82) Antibody. It is currently untested **Product** 

**Comments** in blocking applications.

**Peptide** 

**CVDDSQDYYVGKKN** Sequence

Shipping

Ambient Instructions

**Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring\_medium.png Grade