

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP07069 - IMMUNIZING PEPTIDE

SKU: EBP07069

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage Shipped at ambient temperature, store at -20°C.

Accession

NP_002037.2; NP_001243728.1 ID

The peptide was used in the production of EB07069 - Goat Anti-GAPDH (Internal) Antibody. It is currently untested **Product**

Comments in blocking applications.

Peptide

C-GVNHEKYDNSLK Sequence

Shipping

Ambient Instructions

Product

 $https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png$ Grade