

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP07078 - IMMUNIZING PEPTIDE

SKU: EBP07078

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_005553.2; NP_061486.2 ID

 $\label{eq:lower_lower} \textbf{Immunogen} \ \ \ \ \ \ \ \ \ Peptide with sequence C-QKVDTRAKNAGVTIQD, from the internal region of the protein sequence according to NP_005553.2; NP_061486.2.$

The peptide was used in the production of EB07078 - Goat Anti-LAMC2 Antibody. It is currently untested in **Product**

Comments blocking applications.

Peptide

C-QKVDTRAKNAGVTIQD Sequence

Shipping

Ambient Instructions

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade