

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP07080 - IMMUNIZING PEPTIDE**

**SKU:** EBP07080

## **SPECIFICATIONS**

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$ 

**Unit Size** 100 μα

Shipped at ambient temperature, store at -20°C.

NP\_004947.2; NP\_001156619.1; NP\_001156620.1; NP\_001156621.1; NP\_001156622.1; NP\_001156623.1; Accession

NP\_001156624.1NP\_001357485.1 ID

Peptide with sequence PLKIKKEPHSPCSE, from the internal region of the protein sequence according to

Immunogen NP 004947.2; NP 001156619.1; NP 001156620.1; NP 001156621.1; NP 001156622.1; NP 001156623.1;

NP 001156624.1NP 001357485.1.

The peptide was used in the production of EB07080 - Goat Anti-ETV1 Antibody. It is currently untested in **Product** 

Comments blocking applications.

**Peptide** 

**PLKIKKEPHSPCSE** Sequence

Shipping

Instructions Ambient

**Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png Grade