

EBP07080 - IMMUNIZING PEPTIDE

SKU: EBP07080



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_004947.2; NP_001156619.1; NP_001156620.1; NP_001156621.1; NP_001156622.1; NP_001156623.1; NP_001156624.1NP_001357485.1
Immunogen	Peptide with sequence PLKIKKEPHSPCSE, from the internal region of the protein sequence according to NP_004947.2; NP_001156619.1; NP_001156620.1; NP_001156621.1; NP_001156622.1; NP_001156623.1; NP_001156624.1NP_001357485.1.
Product Comments	The peptide was used in the production of EB07080 - Goat Anti-ETV1 Antibody. It is currently untested in blocking applications.
Peptide Sequence	PLKIKKEPHSPCSE
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)