



Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP07081 - IMMUNIZING PEPTIDE**

**SKU:** EBP07081

## **SPECIFICATIONS**

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

**Unit Size** 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

**Accession** 

 $NP\_891548.1; \ NP\_004440.1; \ NP\_001129626.1; \ NP\_001129627.1; \ NP\_001230358.1 |$ ID

Peptide with sequence C-PNTRLPTSHMPSH, from the C Terminus of the protein sequence according to NP\_891548.1;

Immunogen NP\_004440.1; NP\_001129626.1; NP\_001129627.1; NP\_001230358.1|.

The peptide was used in the production of EB07081 - Goat Anti-ERG Antibody. It is currently untested in blocking **Product** 

Comments applications.

**Peptide** 

C-PNTRLPTSHMPSH Sequence

Shipping

Instructions Ambient

**Product** 

 $https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_plus\_medium.png$ Grade