

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP07082 - IMMUNIZING PEPTIDE**

**SKU:** EBP07082

## **SPECIFICATIONS**

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

**Unit Size** 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

**Accession** 

NP\_031372.2 ID

 $\begin{tabular}{ll} \textbf{Immunogen} \\ \textbf{NP\_031372.2}. \\ \end{tabular} Peptide with sequence C-QDAEVESSAKSGK, from the C Terminus of the protein sequence according to NP\_031372.2. \\ \end{tabular}$ 

**Product** The peptide was used in the production of EB07082 - Goat Anti-OGFR Antibody. It is currently untested in blocking

Comments applications.

**Peptide** 

C-QDAEVESSAKSGK Sequence

**Shipping** 

Ambient Instructions

**Product** 

 $https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_plus\_medium.png$ Grade