

EBP07086 - IMMUNIZING PEPTIDE

SKU: EBP07086



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_004405.3; NP_981962.1; NP_981963.1
Immunogen	Peptide with sequence HVCESDQEKEE, from the internal region of the protein sequence according to NP_004405.3; NP_981962.1; NP_981963.1.
Product Comments	The peptide was used in the production of EBP07086 - Goat Anti-Calciressin 1 / DSCR1 Antibody. It is currently untested in blocking applications.
Peptide Sequence	HVCESDQEKEE
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png