

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP07088 - IMMUNIZING PEPTIDE

SKU: EBP07088

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μα

Storage

Shipped at ambient temperature, store at -20°C. Instructions

Accession

NP 783195.2; NP 006522.2, NP 001340507.1; NP 001340498.1 ID

Peptide with sequence C-KKRIDEDDFADEE, from the C Terminus of the protein sequence according to Immunogen

NP 783195.2; NP 006522.2, NP 001340507.1; NP 001340498.1.

The peptide was used in the production of EB07088 - Goat Anti-IFT88 / Polaris Antibody. It is currently untested **Product**

Comments in blocking applications.

Peptide

C-KKRIDEDDFADEE Sequence

Shipping

Instructions Ambient

Product https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade