

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP07091 - IMMUNIZING PEPTIDE

SKU: EBP07091

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_004515.2; NP_001015002.1; NP_001026973.1 ID

Peptide with sequence C-RVSSRKRHPAGPPGE, from the internal region of the protein sequence according to Immunogen

NP_004515.2; NP_001015002.1; NP_001026973.1.

The peptide was used in the production of EB07091 - Goat Anti-LLGL2 Antibody. It is currently untested in blocking Product

Comments applications.

Peptide

C-RVSSRKRHPAGPPGE Sequence

Ambient Instructions Shipping

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade