

EBP07111 - IMMUNIZING PEPTIDE

SKU: EBP07111



SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 µg

Storage Instructions Shipped at ambient temperature, store at -20°C.

Accession ID NP_000556.1

Immunogen Peptide with sequence C-DNTSSHNRPDARD, from the internal region of the protein sequence according to NP_000556.1.

Product Comments The peptide was used in the production of EBP07111 - Goat Anti-IL6R / CD126 (isoform 1) Antibody. It is currently untested in blocking applications.

Peptide Sequence C-DNTSSHNRPDARD

Shipping Instructions Ambient

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

DOCUMENTS

- [Data Sheet](#)