

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP07116 - IMMUNIZING PEPTIDE

SKU: EBP07116

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession NP_001012975.1; NP_001012976.1; NP_001012977.1; NP_061132.2; NP_001158009.1; NP_001158010.1;

NP_001158011.1; NP_001158012.1 and more ID

Peptide with sequence C-GEESHHSESRARQ, from the internal region of the protein sequence according to

Immunogen NP_001012975.1; NP_001012976.1; NP_001012977.1; NP_061132.2; NP_001158009.1; NP_001158010.1;

NP_001158011.1; NP_001158012.1 and more.

The peptide was used in the production of EB07116 - Goat Anti-DISC1 Antibody. It is currently untested in blocking **Product**

Comments applications.

Peptide

C-GEESHHSESRARQ

Sequence

Shipping Ambient Instructions

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade