

EBP07116 - IMMUNIZING PEPTIDE

SKU: EBP07116



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_001012975.1; NP_001012976.1; NP_001012977.1; NP_061132.2; NP_001158009.1; NP_001158010.1; NP_001158011.1; NP_001158012.1 and more
Immunogen	Peptide with sequence C-GEESHSESRRARQ, from the internal region of the protein sequence according to NP_001012975.1; NP_001012976.1; NP_001012977.1; NP_061132.2; NP_001158009.1; NP_001158010.1; NP_001158011.1; NP_001158012.1 and more.
Product Comments	The peptide was used in the production of EBP07116 - Goat Anti-DISC1 Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-GEESHSESRRARQ
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png