

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP07119 - IMMUNIZING PEPTIDE

SKU: EBP07119

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_055512.1 ID

Immunogen Peptide with sequence CNVFKSKKKSQ, from the internal region of the protein sequence according to NP_055512.1.

Product The peptide was used in the production of EB07119 - Goat Anti-CAPON / NOS1AP (internal) Antibody. It is currently

Comments untested in blocking applications.

Peptide

Sequence

CNVFKSKKKSQ

Shipping Ambient Instructions

Product

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png