

EBP07119 - IMMUNIZING PEPTIDE

SKU: EBP07119



SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 µg

Storage Instructions Shipped at ambient temperature, store at -20°C.

Accession ID NP_055512.1

Immunogen Peptide with sequence CNVFKSKKKSQ, from the internal region of the protein sequence according to NP_055512.1.

Product Comments The peptide was used in the production of EBP07119 - Goat Anti-CAPON / NOS1AP (internal) Antibody. It is currently untested in blocking applications.

Peptide Sequence CNVFKSKKKSQ

Shipping Instructions Ambient

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png