

EBP07152 - IMMUNIZING PEPTIDE

SKU: EBP07152



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_000948.2; NP_942005.1; NP_942006.1; NP_942007.1; NP_942012.1; NP_942013.1; NP_942008.1; NP_942009.1; NP_942010.1; NP_942011.1
Immunogen	Peptide with sequence HCKTHTEKQKECN, from the internal region of the protein sequence according to NP_000948.2; NP_942005.1; NP_942006.1; NP_942007.1; NP_942012.1; NP_942013.1; NP_942008.1; NP_942009.1; NP_942010.1; NP_942011.1.
Product Comments	The peptide was used in the production of EB07152 - Goat Anti-Prostaglandin E Receptor 3 Antibody. It is currently untested in blocking applications.
Peptide Sequence	HCKTHTEKQKECN
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

DOCUMENTS

- [Data Sheet](#)