



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP07171 - IMMUNIZING PEPTIDE

SKU: EBP07171

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage Shipped at ambient temperature, store at -20°C.

Accession

NP_066961.2 ID

Immunogen Peptide with sequence C-KMEAKNSYVRDD, from the C Terminus of the protein sequence according to NP_066961.2. **Product** The peptide was used in the production of EB07171 - Goat Anti-TRAF2 Antibody. It is currently untested in blocking

Comments applications.

Peptide

C-KMEAKNSYVRDD Sequence

Shipping

Ambient Instructions

Product

 $https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png$ Grade