

# EBP07177 - IMMUNIZING PEPTIDE

SKU: EBP07177



---

## SPECIFICATIONS

**Formulation** 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

**Unit Size** 100 µg

**Storage Instructions** Shipped at ambient temperature, store at -20°C.

**Accession ID** NP\_006425.2; NP\_056200.1; NP\_001030126.1; NP\_001030127.1; NP\_001030128.1; NP\_079267.1; NP\_001030129.1; NP\_001277227.2

**Immunogen** Peptide with sequence C-SPKSEDDDDSDLYSP, from the internal region of the protein sequence according to NP\_006425.2; NP\_056200.1; NP\_001030126.1; NP\_001030127.1; NP\_001030128.1; NP\_079267.1; NP\_001030129.1; NP\_001277227.2.

**Product Comments** The peptide was used in the production of EB07177 - Goat Anti-SORBS1 / ponsin (internal) Antibody. It is currently untested in blocking applications.

**Peptide Sequence** C-SPKSEDDDDSDLYSP

**Shipping Instructions** Ambient

**Product Grade** [https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\\_medium.png](https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png)

## DOCUMENTS

- [Data Sheet](#)