

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP07177 - IMMUNIZING PEPTIDE

SKU: EBP07177

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μα

Accession

Shipped at ambient temperature, store at -20°C.

NP_006425.2; NP_056200.1; NP_001030126.1; NP_001030127.1; NP_001030128.1; NP_079267.1; ID

NP_001030129.1; NP_001277227.2

Peptide with sequence C-SPKSEDDDSDLYSP, from the internal region of the protein sequence according to Immunogen NP 006425.2; NP 056200.1; NP 001030126.1; NP 001030127.1; NP 001030128.1; NP 079267.1;

NP 001030129.1; NP 001277227.2.

The peptide was used in the production of EB07177 - Goat Anti-SORBS1 / ponsin (internal) Antibody. It is **Product**

Comments currently untested in blocking applications.

Peptide

C-SPKSEDDDSDLYSP Sequence

Shipping

Instructions Ambient

Product https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade