

EBP07181 - IMMUNIZING PEPTIDE

SKU: EBP07181



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_004161.4
Immunogen	Peptide with sequence C-DNEDSDTKKSYVN, from the internal region of the protein sequence according to NP_004161.4.
Product	The peptide was used in the production of EB07181 - Goat Anti-EAAC1/ SLC1A1 / EAAT3 Antibody. It is
Comments	currently untested in blocking applications.
Peptide Sequence	C-DNEDSDTKKSYVN
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png