

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP07196 - IMMUNIZING PEPTIDE

SKU: EBP07196

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Shipped at ambient temperature, store at -20°C.

Accession

 $NP_057353.1; \ NP_001124185.1; \ NP_001124186.1; \ NP_001159591.1$ ID

Peptide with sequence C-QHEQRKEQEPKRPH, from the internal region of the protein sequence according to Immunogen NP_057353.1; NP_001124185.1; NP_001124186.1; NP_001159591.1.

The peptide was used in the production of EB07196 - Goat Anti-LEF1 Antibody. It is currently untested in blocking **Product**

Comments applications.

Peptide

C-QHEQRKEQEPKRPH Sequence

Shipping

Ambient Instructions

Product

 $https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png$ Grade