



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP07207 - IMMUNIZING PEPTIDE

SKU: EBP07207

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage Instructions Shipped at ambient temperature, store at -20°C.

Accession NP_006658.1 ID

Immunogen Peptide with sequence C-EPKDESARKND, from the C Terminus of the protein sequence according to NP_006658.1.

Product The peptide was used in the production of EB07207 - Goat Anti-PGRMC1 / MPR Antibody. It is currently untested in

Comments blocking applications.

Peptide

C-EPKDESARKND Sequence

Shipping

Ambient Instructions

Product

 $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. amazonaws. com/assets/products/image/elite_plus_medium.png$ Grade