

EBP07209 - IMMUNIZING PEPTIDE

SKU: EBP07209



SPECIFICATIONS

| | |
|------------------------------|---|
| Formulation | 100µg of Lyophilised peptide, formulation De-ionized (DI) water. |
| Unit Size | 100 µg |
| Storage | |
| Instructions | Shipped at ambient temperature, store at -20°C. |
| Accession ID | NP_066964.1 |
| Immunogen | Peptide with sequence C-SLAKKDEKDTLED, from the internal region of the protein sequence according to NP_066964.1. |
| Product Comments | The peptide was used in the production of EB07209 - Goat Anti-Ku80 / XRCC5 Antibody. It is currently untested in blocking applications. |
| Peptide Sequence | C-SLAKKDEKDTLED |
| Shipping Instructions | Ambient |
| Product Grade | https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png |

DOCUMENTS

- [Data Sheet](#)