

# EBP07220 - IMMUNIZING PEPTIDE

**SKU:** EBP07220



## SPECIFICATIONS

<b>Formulation</b>	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
<b>Unit Size</b>	100 µg
<b>Storage Instructions</b>	Shipped at ambient temperature, store at -20°C.
<b>Accession ID</b>	NP_003409.1;; NP_001120664.1; NP_001120665.1;NP_001120666.1
<b>Immunogen</b>	Peptide with sequence GESGHLARECTIE, from the C Terminus of the protein sequence according to NP_003409.1;; NP_001120664.1; NP_001120665.1;NP_001120666.1.
<b>Product Comments</b>	The peptide was used in the production of EB07220 - Goat Anti-ZNF9 / CNBP Antibody. It is currently untested in blocking applications.
<b>Peptide Sequence</b>	GESGHLARECTIE
<b>Shipping Instructions</b>	Ambient
<b>Product Grade</b>	<a href="https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png">https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png</a>

## DOCUMENTS

- [Data Sheet](#)