

EBP07233 - IMMUNIZING PEPTIDE

SKU: EBP07233



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_066124.1; NP_065681.1; NP_001342145.1
Immunogen	Peptide with sequence NDTKALRRPKCAELH, from the internal region of the protein sequence according to NP_066124.1; NP_065681.1; NP_001342145.1.
Product	The peptide was used in the production of EBP07233 - Goat Anti-RET Antibody. It is currently untested in
Comments	blocking applications.
Peptide Sequence	NDTKALRRPKCAELH
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png