

EBP07239 - IMMUNIZING PEPTIDE

SKU: EBP07239



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	Shipped at ambient temperature, store at -20°C.
Instructions	
Accession ID	NP_000638.1; NP_000639.1
Immunogen	Peptide with sequence RFIRNTNESGEEVT-C, from the N Terminus of the protein sequence according to NP_000638.1; NP_000639.1.
Product	The peptide was used in the production of EB07239 - Goat Anti-CCR2 Antibody. It is currently untested in blocking applications.
Comments	
Peptide Sequence	RFIRNTNESGEEVT-C
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png