

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP07241 - IMMUNIZING PEPTIDE

SKU: EBP07241

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μα

Shipped at ambient temperature, store at -20°C.

Accession

NP_001361195.1 ID

Peptide with sequence CDLNCRIHKKSRNK, from the internal region of the protein sequence according to

Immunogen NP 005028.5; NP 056953.2; NP 001341597.1; NP 001361194.1; NP 001341599.1; NP 001317544.2;

NP 001361195.1.

The peptide was used in the production of EB07241 - Goat Anti-PPARG Antibody. It is currently untested in **Product**

Comments blocking applications.

Peptide

CDLNCRIHKKSRNK Sequence

Shipping

Instructions Ambient

Product https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade