

# EBP07241 - IMMUNIZING PEPTIDE

**SKU:** EBP07241



---

## SPECIFICATIONS

**Formulation** 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

**Unit Size** 100 µg

**Storage Instructions** Shipped at ambient temperature, store at -20°C.

**Accession ID** NP\_005028.5; NP\_056953.2; NP\_001341597.1; NP\_001361194.1; NP\_001341599.1; NP\_001317544.2; NP\_001361195.1

**Immunogen** Peptide with sequence CDLNCRIHKKSRNK, from the internal region of the protein sequence according to NP\_005028.5; NP\_056953.2; NP\_001341597.1; NP\_001361194.1; NP\_001341599.1; NP\_001317544.2; NP\_001361195.1.

**Product Comments** The peptide was used in the production of EB07241 - Goat Anti-PPARG Antibody. It is currently untested in blocking applications.

**Peptide Sequence** CDLNCRIHKKSRNK

**Shipping Instructions** Ambient

**Product Grade** [https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\\_medium.png](https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png)