

EBP07254 - IMMUNIZING PEPTIDE

SKU: EBP07254



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_570722.2; NP_001028731.1; NP_001028732.1
Immunogen	Peptide with sequence PLQVSDYSWQTKT-C, from the N Terminus of the protein sequence according to NP_570722.2; NP_001028731.1; NP_001028732.1.
Product Comments	The peptide was used in the production of EB07254 - Goat Anti-DYX1C1 / EKN1 Antibody. It is currently untested in blocking applications.
Peptide Sequence	PLQVSDYSWQTKT-C
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)