

EBP07291 - IMMUNIZING PEPTIDE

SKU: EBP07291



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_115498.2; NP_898861.1
Immunogen	Peptide with sequence C-KTLSDKSREAKVK, from the internal region of the protein sequence according to NP_115498.2; NP_898861.1.
Product	The peptide was used in the production of EBP07291 - Goat Anti-DTNBP1 / Dysbindin Antibody. It is currently
Comments	untested in blocking applications.
Peptide Sequence	C-KTLSDKSREAKVK
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png