

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP07322 - IMMUNIZING PEPTIDE**

**SKU:** EBP07322

## **SPECIFICATIONS**

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

**Unit Size** 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

**Accession** 

NP\_055395.2 ID

Immunogen Peptide with sequence C-KGKTLEQITAH, from the C Terminus of the protein sequence according to NP\_055395.2.

Product The peptide was used in the production of EB07322 - Goat Anti-GLUT8 / SLC2A8 Antibody. It is currently untested

**Comments** in blocking applications.

**Peptide** 

Sequence

C-KGKTLEQITAH

**Shipping** Ambient Instructions

**Product** 

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring\_medium.png