

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP07328 - IMMUNIZING PEPTIDE

SKU: EBP07328

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Shipped at ambient temperature, store at -20°C.

Accession

ID

NP_998786.3 NP_000860.3; NP_001155244.1

The peptide was used in the production of EB07328 - Goat Anti-Serotonin receptor 3A / HTR3A Antibody. It is **Product**

Comments currently untested in blocking applications.

Peptide

Sequence Shipping

C-GPQDFEKSPRDR

Ambient Instructions

Product

 $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. amazonaws. com/assets/products/image/elite_plus_medium.png$ Grade