

EBP07333-B - IMMUNIZING PEPTIDE

SKU: EBP07333-B



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	Shipped at ambient temperature, store at -20°C.
Instructions	
Accession ID	NP_000542.1; NP_937799.1; XP_011532380.1
Immunogen	Peptide with sequence C-RSLVKPENYRRLD., from the internal region of the protein sequence according to NP_000542.1; NP_937799.1; XP_011532380.1.
Product	The peptide was used in the production of EB07333-B - Goat Anti-VHL, Biotinylated Antibody. It is currently untested in blocking applications.
Comments	
Peptide Sequence	C-RSLVKPENYRRLD.
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)