

EBP07341 - IMMUNIZING PEPTIDE

SKU: EBP07341



SPECIFICATIONS

| | |
|-------------------------|---|
| Formulation | 100µg of Lyophilised peptide, formulation De-ionized (DI) water. |
| Unit Size | 100 µg |
| Storage | |
| Instructions | Shipped at ambient temperature, store at -20°C. |
| Accession ID | NP_000823.4; NP_067544.1; NP_015566.1 |
| Immunogen | Peptide with sequence C-TQERVNNSNKKE, from the internal region of the protein sequence according to NP_000823.4; NP_067544.1; NP_015566.1. |
| Product | The peptide was used in the production of EB07341 - Goat Anti-NMDA receptor 1 / GRIN1 Antibody. It is currently untested in blocking applications. |
| Comments | |
| Peptide Sequence | C-TQERVNNSNKKE |
| Shipping | |
| Instructions | Ambient |
| Product Grade | https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png |

DOCUMENTS

- [Data Sheet](#)