

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP07378 - IMMUNIZING PEPTIDE

SKU: EBP07378

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Shipped at ambient temperature, store at -20°C.

Accession

NP_003097.1 ID

 $\label{eq:lower_lower_lower} \textbf{Immunogen} \ \ \ \ \ \ \ \ Peptide \ with \ sequence \ YNMMETELKPPGPQQ-C, from \ the \ N \ Terminus \ of the protein \ sequence \ according to \ NP_003097.1.$

Product The peptide was used in the production of EB07378 - Goat Anti-SOX2 Antibody. It is currently untested in blocking

Comments applications.

Peptide

YNMMETELKPPGPQQ-C Sequence

Shipping

Ambient Instructions

Product

 $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. amazonaws. com/assets/products/image/elite_plus_medium.png$ Grade