

# EBP07379 - IMMUNIZING PEPTIDE

**SKU:** EBP07379



## SPECIFICATIONS

|                              |                                                                                                                                                                                                                                             |
|------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Formulation</b>           | 100µg of Lyophilised peptide, formulation De-ionized (DI) water.                                                                                                                                                                            |
| <b>Unit Size</b>             | 100 µg                                                                                                                                                                                                                                      |
| <b>Storage Instructions</b>  | Shipped at ambient temperature, store at -20°C.                                                                                                                                                                                             |
| <b>Accession ID</b>          | NP_644805.1; NP_003141.2; NP_998827.1                                                                                                                                                                                                       |
| <b>Immunogen</b>             | Peptide with sequence CRPESQEHPEADP, from the internal region of the protein sequence according to NP_644805.1; NP_003141.2; NP_998827.1.                                                                                                   |
| <b>Product Comments</b>      | The peptide was used in the production of EB07379 - Goat Anti-STAT3 Antibody. It is currently untested in blocking applications.                                                                                                            |
| <b>Peptide Sequence</b>      | CRPESQEHPEADP                                                                                                                                                                                                                               |
| <b>Shipping Instructions</b> | Ambient                                                                                                                                                                                                                                     |
| <b>Product Grade</b>         | <a href="https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png">https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png</a> |

## DOCUMENTS

- [Data Sheet](#)