

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP07382 - IMMUNIZING PEPTIDE

SKU: EBP07382

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

100 μg **Unit Size**

Storage

Shipped at ambient temperature, store at -20°C. Instructions

Accession

NP 443180.2; NP 001353107.1 ID

Peptide with sequence ARSTDSLDGPGE-C, from the N Terminus of the protein sequence according to Immunogen NP_443180.2; NP_001353107.1.

The peptide was used in the production of EB07382 - Goat Anti-ARHGAP33 / TCGAP (N Terminus) Antibody. It is **Product**

Comments currently untested in blocking applications.

Peptide ARSTDSLDGPGE-C

Sequence

Shipping Instructions Ambient

Product https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade