

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP07386 - IMMUNIZING PEPTIDE

SKU: EBP07386

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_003099.1 ID

 $\begin{tabular}{ll} \textbf{Immunogen} \\ \textbf{NP_003099.1}. \\ \end{tabular} \begin{tabular}{ll} \textbf{Peptide with sequence KTSKGSSKKCGK, from the internal region of the protein sequence according to NP_003099.1. \\ \end{tabular}$

The peptide was used in the production of EB07386 - Goat Anti-SOX11 / SRY-box11 Antibody. It is currently **Product**

Comments untested in blocking applications.

Peptide

KTSKGSSKKCGK Sequence

Shipping

Ambient Instructions

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade