

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP07387 - IMMUNIZING PEPTIDE

SKU: EBP07387

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage Instructions Shipped at ambient temperature, store at -20°C.

Accession NP_000705.2 ID

The peptide was used in the production of EB07387 - Goat Anti-TSPO / PBR (isoform) Antibody. It is currently **Product**

Comments untested in blocking applications.

Peptide

C-RDNHGWRGGRRLPE Sequence

Shipping

Ambient Instructions

Product

 $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. amazonaws. com/assets/products/image/elite_plus_medium.png$ Grade