

EBP07420 - IMMUNIZING PEPTIDE

SKU: EBP07420



SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 µg

Storage Instructions Shipped at ambient temperature, store at -20°C.

Accession ID NP_000785.1

Immunogen Peptide with sequence C-DYDTDVSLEKIQ, from the C Terminus of the protein sequence according to NP_000785.1.

Product Comments The peptide was used in the production of EBP07420 - Goat Anti-Dopamine receptor D1 / DRD1 Antibody. It is currently untested in blocking applications.

Peptide Sequence C-DYDTDVSLEKIQ

Shipping Instructions Ambient

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

DOCUMENTS

- [Data Sheet](#)