

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP07446 - IMMUNIZING PEPTIDE

SKU: EBP07446

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μα

Storage

Shipped at ambient temperature, store at -20°C. Instructions

Accession

NP_000676.1; NP_004826.2; NP_033611.1; NP_114038.1; NP_114438.1 ID

Peptide with sequence C-EIQKNKPRNDDIFK, from the internal region of the protein sequence according to Immunogen

NP 000676.1; NP 004826.2; NP 033611.1; NP 114038.1; NP 114438.1.

The peptide was used in the production of EB07446 - Goat Anti-AGTR1 / AT1 Antibody. It is currently untested **Product**

Comments in blocking applications.

Peptide

C-EIQKNKPRNDDIFK Sequence

Shipping

Instructions Ambient

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png