

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP07454 - IMMUNIZING PEPTIDE

SKU: EBP07454

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μα

Shipped at ambient temperature, store at -20°C.

NP_663731.1; NP_001248356.1; NP_001248357.1; NP_001336129.1; NP_001336131.1; NP_001336133.1; Accession

NP_001336136.1; NP_001336137.1 ID

Peptide with sequence C-NKTDSPSRKRDK, from the internal region of the protein sequence according to

Immunogen NP 663731.1; NP 001248356.1; NP 001248357.1; NP 001336129.1; NP 001336131.1; NP 001336133.1;

NP 001336136.1; NP 001336137.1.

The peptide was used in the production of EB07454 - Goat Anti-Lipin 1 (Internal) Antibody. It is currently **Product**

Comments untested in blocking applications.

Peptide

C-NKTDSPSRKRDK Sequence

Shipping

Instructions Ambient

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade