

EBP07470 - IMMUNIZING PEPTIDE

SKU: EBP07470



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_005616.2; NP_001007068.1; NP_001007069; NP_001007070.1; NP_001007071.1
Immunogen	Peptide with sequence SLEDLKVDKVIQAQ-C, from the N Terminus of the protein sequence according to NP_005616.2; NP_001007068.1; NP_001007069; NP_001007070.1; NP_001007071.1.
Product Comments	The peptide was used in the production of EB07470 - Goat Anti-Syntenin / SDCBP Antibody. It is currently untested in blocking applications.
Peptide Sequence	SLEDLKVDKVIQAQ-C
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)