

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP07483 - IMMUNIZING PEPTIDE

SKU: EBP07483

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_001845.3; NP_542196.2; NP_542197.2 ID

 $\label{eq:local_$

The peptide was used in the production of EB07483 - Goat Anti-COL11A1 Antibody. It is currently untested in Product

Comments blocking applications. **Peptide**

C-EDYDSQRKNSEDTLY Sequence

Shipping

Ambient Instructions

Product https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade