

EBP07489 - IMMUNIZING PEPTIDE

SKU: EBP07489



SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 µg

Storage Instructions Shipped at ambient temperature, store at -20°C.

Accession ID NP_055081.1; NP_001007470.1

Immunogen Peptide with sequence C-RSPAQEDGKVYIN, from the C Terminus of the protein sequence according to NP_055081.1; NP_001007470.1.

Product Comments The peptide was used in the production of EBP07489 - Goat Anti-DAP10 / HCST Antibody. It is currently untested in blocking applications.

Peptide Sequence C-RSPAQEDGKVYIN

Shipping Instructions Ambient

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)