

EBP07498 - IMMUNIZING PEPTIDE

SKU: EBP07498



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_004903.2; NP_919436.1; NP_919437.1; NP_919438.1; NP_001013424.1
Immunogen	Peptide with sequence C-YGNYESKKHKQG, from the internal region of the protein sequence according to NP_004903.2; NP_919436.1; NP_919437.1; NP_919438.1; NP_001013424.1.
Product Comments	The peptide was used in the production of EB07498 - Goat Anti-KRIT1 / CCM1 Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-YGNYESKKHKQG
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

DOCUMENTS

- [Data Sheet](#)