

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP07498 - IMMUNIZING PEPTIDE

SKU: EBP07498

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_004903.2; NP_919436.1; NP_919437.1; NP_919438.1; NP_001013424.1 ID

The peptide was used in the production of EB07498 - Goat Anti-KRIT1 / CCM1 Antibody. It is currently untested in Product

Comments blocking applications.

Peptide

C-YGNYESKKHKQG Sequence

Shipping

Ambient Instructions

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade