

EBP07500 - IMMUNIZING PEPTIDE

SKU: EBP07500



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_001448.2; NP_001157789.1; NP_001157790.1; NP_001157791.1
Immunogen	Peptide with sequence KITDDSRRC SQVK, from the internal region of the protein sequence according to NP_001448.2; NP_001157789.1; NP_001157790.1; NP_001157791.1.
Product Comments	The peptide was used in the production of EBP07500 - Goat Anti-FLNB Antibody. It is currently untested in blocking applications.
Peptide Sequence	KITDDSRRC SQVK
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

DOCUMENTS

- [Data Sheet](#)