

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP07500 - IMMUNIZING PEPTIDE

SKU: EBP07500

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_001448.2; NP_001157789.1; NP_001157790.1; NP_001157791.1 ID

Peptide with sequence KITDDSRRCSQVK, from the internal region of the protein sequence according to Immunogen

 $NP_001448.2; \ NP_001157789.1; \ NP_001157790.1; \ NP_001157791.1.$

The peptide was used in the production of EB07500 - Goat Anti-FLNB Antibody. It is currently untested in blocking Product

Comments applications.

Peptide

KITDDSRRCSQVK Sequence

Instructions Ambient Shipping

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade